

TERUMO PENPOL REVIEW

ON BLOOD MANAGEMENT SYSTEMS

JAN-MAR **17** NUMBER **89**

Copyright :

© Asian Journal of Transfusion Science

This is an openaccess article distributed under the terms of the Creative Commons Attribution Non commercial Share Alike 3.0 Unported, which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

Asian J Transfus Sci. 2014 JulDec; 8(2): 89–91.

For Private Circulation For Educational Purpose Only

Published in every quarter by TERUMO PENPOL Private Limited

TERUMO PENPOL®

Experience of buffy coat pooling of platelets as a supportive care in thrombocytopenic dengue patients : A prospective study

Kabita Chatterjee, Poonam Coshic, Mayuri Borgohain, and Nitin Agarwal

Department of Transfusion Medicine, Blood Bank (Main), All India Institute of Medical Sciences, Ansari Nagar, New Delhi, India

Correspondence to: Dr. Kabita Chatterjee, Room No. 1, Blood Bank (Main), All India Institute of Medical Sciences, Ansari Nagar, New Delhi – 110 029, India. Email : kabita.chatterjee@gmail.com

Abstract

Random donor platelet (RDP) is not sufficient to improve the platelet count in most thrombocytopenic patients. Single donor platelet (SDP) or buffy coat pooled platelet (BCPP) are the two choices to provide a full therapeutic dose of platelets. However, there are constraints in the preparation of SDP due to stringent donor selection procedure, time required for

procedure, and need of special expensive equipments and kits. BCPP is widely practiced, especially in the European countries, since 1995. In India, we decided to adopt the procedure of buffy coat pooling of platelets, especially for economically backward patients and for emergencies. This study was prospectively conducted from September 2009 to September 2010. A total of 129 units of BCPP [tested

prior for viral markers by enzyme linked immunosorbent assay (ELISA) and individual donor nucleic acid amplification test (IDNAT)] were issued to 129 patients suffering from dengue and were included in this study. For comparison between efficacy of SDP and BCCP, patients were divided into two groups of 50 each. The posttransfusion platelet counts of the patients were noted after 2 hours of transfusion for each type of component. The platelet yield varied from 2.5 to 4.4 $\times 10^{10}$ in BCCP samples. The samples analyzed were sterile without any contamination. The different biochemical parameters were analyzed in detail. The observed posttransfusion platelet recovery and corrected count increment (CCI) at 1 hour and 24 hours after BCCP transfusion were similar to that after SDP transfusion. Hence, we concluded that BCCP can be a low cost alternative to SDP in the times of emergencies like dengue and nonaffordability by the patient for SDP.

Keywords: buffy coat, buffy coat pooling of platelets, random donor platelets, single donor platelets

Introduction

Therapeutic or prophylactic use of platelet concentrates is vital for patients with thrombocytopenia due to dengue fever, intensive chemotherapy, and other illnesses. Dengue is one of the major diseases in and around Delhi. Dengue infections are seen every year, thus making it an endemic disease.[1] The mechanisms underlying the bleeding in dengue are multiple. These are vasculopathy, thrombopathy, and disseminated intravascular coagulopathy. Thrombopathy consists of thrombocytopenia and platelet dysfunction.[2] Bleeding occurs more often in patients with severe thrombocytopenia. Highrisk patients having platelet count $<20,000/\text{mm}^3$ and risk of bleeding require urgent platelet transfusion.[3]

Platelet concentrates can be obtained by platelet rich plasma (PRP) method, buffy coat method, and apheresis procedure using cell separator machines.[4,5,6] New trends in the preparation of platelets showed a shift toward the use of the buffy coat method. It has been suggested that this method causes less platelet

activation and damage during platelet preparation. Platelet preparation by buffy coat pooling method is practiced in European countries, US, Latin America, and in some Asian countries.[4]

The buffy coat pooling method of harvesting platelets, which is equivalent to single donor platelet (SDP) apheresis, helps to improve platelet yield, meets emergency requirement for platelet, reduces cost to poor patients, and reduces leukocyte contamination in platelet.[4] In the present study, 129 pooled platelet units were prepared and issued to 129 dengue patients, and a study was conducted on platelet yield, storage parameters, and also on sterility parameters on all issued units. Observational study was also conducted on transfused patients and the posttransfusion platelet increment was compared with that of apheresis platelets.

Material and Methods

Buffy coat pooled platelet (BCPP) preparation

The four buffy coat bags prepared by centrifugation of donated wholeblood were hung vertically from a stand for 4 hours at 22°C until the

confirmed reports of ABO blood grouping and transfusion transmitted infections (TTI) by enzyme linked immunosorbent assay (ELISA) and by individual donor nucleic acid amplification testing (IDNAT) were obtained. The buffy coat bags were connected serially using sterile connecting device(TSCD II, Terumo medical corp, NJ). One plasma unit (200 ml) of the same blood group was connected to the topmost bag. Plasma going down due to gravity through each bag (which were connected to each other), rinsed each bag in the process, and got collected in the bottom most bag. The bottommost bag containing 400 ml of fluid (plasma + 4 buffy coat bags) was centrifuged (Cryofuge 6000i, Heraeus, Thermo Scientific) at 1240 rpm for 10 min. After centrifugation, platelets were separated from platelet pool. The Log 3 leucoreduced platelets were collected in a special type 1000 ml platelet bag and stored at 22°C ± 2°C in platelet agitator/incubator for 5 days; however, BCPP were issued within 2448 hours of preparation in our study.

Sampling

A total of 10 ml of platelets were collected into sample bag attached to the main platelet bag and separated maintaining the closed system. This sample bag was stored in similar conditions as the main bag and used for assessment of platelet quality. A total of 129 BCPP samples were analyzed on the first, third, and the fifth days of storage for platelet count (using Beckman Coulter Cell counter), glucose, lactate levels, and pH using blood gas analyzer (Nova). BactAlert

comparison between efficacy of SDP and BCCP, two groups were studied (50 patients who were transfused with SDP and first 50 patients with BCPP). The posttransfusion platelet counts of the patients were noted after 2 hours of transfusion for each type of the component. The CCI[7] and posttransfusion platelet recovery (PPR)[7] were calculated as follows:

Total blood volume was calculated as 75 ml/kg body weight. The Drugs Controller General of India (DCGI) was intimated about

$$CCI = \frac{(\text{Post transfusion count} - \text{pre-transfusion count}) \times BSA (m^2)}{\text{Number of platelets transfused}}$$

$$PPR = \frac{(\text{Post transfusion count} - \text{pre transfusion count}) \times \text{Total blood volume}}{\text{Number of platelets transfused}}$$

(Biomerieux) was also used to check bacterial contamination. Residual platelet counts after filtration were done by Coulter cell counter, and leucocyte count was done by Coulter cell counter as well as Nageotte chamber.

SDPs were collected using Hemonetics MCS+ and Fresenius COM.TEC machines on healthy donors with hemoglobin >12.5% and weight >50 kg. For

the process of buffy coat pooling of platelets, and we were granted the license for the pooling of platelets by the Drug Controller authorities.

Results

In the present study, most of the patients belonged to the 20-29 years age group, and the median age was 27 years; 93 patients suffering from dengue who were transfused with BC pooled platelets were males, while 36 were females.

Table 1 shows the platelet counts and other biochemical parameters of BCPP platelets from days 15 of storage. The platelet count did not change much from day 1 to day 5. The viability and metabolic function of platelets were maintained during the 5 days of storage as depicted by the pH and glucose and lactate levels during the storage.

As seen from Table 2, there was not much variation in the platelet count of BCPP and apheresis platelets. The same holds true for the pH values also. By using platelet filter, the leucocyte counts were reduced to Log 3 (1000 times less than the prefiltration value, 0.077×10^6 to 0.001×10^6), but the platelet count changed very slightly. Also, there was no bacterial contamination found in both SDP and BCPP samples by BactAlert done on days 1, 3, 5, and 7 of collection.

Table 3 shows clearly that, in terms of clinical benefit to the patients, BCPP and SDP are not much different. Average increment in platelet count, CCI, and PPR after BCPP transfusion were comparable to that after SDP transfusion.

Discussion

Usage of various platelet products varies greatly between countries and individual institutions.[4] Singh et al.,[5] concluded that, in developing countries, apheresis platelets because of their high cost and more technical expertise required, may be recommended only in selected patients either when PRP platelet concentrates and BCPPs in adequate doses are not available in the inventory mostly due to time constraints.

The results of this study provide evidence that the viability and metabolic function of platelets were maintained as depicted by the glucose and lactate levels during the 5day storage period. Also, the pH range of BC pooled platelets complies with the Council of European Guidelines quality requirements for apheresis platelets (6.47.4).

The average increment in platelet count 2 hours after transfusion of one unit of BCPP is almost equal to that with the transfusion of one unit of SDP. Although, BCPP transfusion involves exposure to more donors than a SDP product, but all the products were found to be

sterile by BactAlert and, thus, were comparable to SDP in these terms. Also, BCPP was prepared only after testing of donors for viral infections by both ELISA and IDNAT, making it further safe for the recipients. Apheresis procedure requires a donor who meets the criteria required for the procedure and also a special kit that cannot be afforded by those with poor financial backgrounds. Most of the patients at our institution needing platelets are below the poverty line (BPL), and treatment affordability is a major concern. As BCPP preparation does not require any such kit, it thereby reduces the cost of platelet therapy.

The preprocedure investigations for SDP need to be carried out (for the SDP donor) prior to the procedure. Our long experience of preparing SDP since 1996 reveals that most of the admitted dengue patients are hemorrhagic type with platelets counts <10,000. In these cases, clinicians demand immediate platelet transfusion and it becomes difficult to wait for 7 to 8 hours for the SDP donor to

be screened and the procedure to be carried out. It is also not easy to generate adequate voluntary apheresis donors to meet the ever increasing demand for SDP due to reasons like time constraints. Our rigorous efforts to create a pool of voluntary apheresis donors has failed to gain momentum in spite of our sincere efforts, and this forced us to explore the possibility of buffy coat pooling of platelets. In order to improve the daytoday patient care services, two pooled buffy coat platelets of each group are prepared daily for emergency purposes at our centre.

Keeping in mind the above constraints and limitations and with the aim to serve the needy patients, we took the initiative to carry out the buffy coat pooling of platelets. We have incorporated certain protocols in our institute to provide better services regarding the issue of platelets which is, in dengue crisis, if the demand for SDP comes in the night hours or if the SDP donor is not available, BCPP is issued.

Conclusions

We found BCPP very valuable during the dengue fever season in 2009-2010.

BCPP can be a good alternative for SDP for correcting thrombocytopenia in cases of emergencies like dengue, nonavailability of donors, and nonaffordability of patients for apheresis.

Footnotes

Source of Support: Nil

Conflict of Interest: None declared

References

- Gupta E, Dar L, Kapoor G, Broor S. The changing epidemiology of dengue in Delhi, India. *Virology*. 2006;3:92. [PMCID: PMC1636631] [PubMed: 17083743]
- Srichaikul T, Nimmannitya S. Haematology in dengue and dengue haemorrhagic fever. *Baillieres Best Pract Res Clin Haematol*. 2000;13:261–76.
- Makroo RN, Raina V, Kumar P, Kanth RK. Role of platelet transfusion in the management of dengue patients in a tertiary care hospital. *Asian J Transfus Sci*. 2007;1:4–7. [PMCID: PMC3168133] [PubMed: 21938225]
- Schrezenmeier H, Seifried E. Buffycoatderived pooled platelet concentrates and apheresis platelet concentrates: Which product type should be preferred? *Vox Sang*. 2010;99:1–15. [PubMed: 20059760]
- Singh RP, Marwaha N, Malhotra P, Dash S. Quality assessment of platelet concentrates prepared by platelet rich plasma platelet concentrate, buffy coat poor platelet concentrate (BCPC) and apheresis PC methods. *Asian J Transfus Sci*. 2009;3:86–94. [PMCID: PMC2920479] [PubMed: 20808653]
- Ralph RV, Scott M. A critical comparison of platelet preparation methods. *Curr Opin Hematol*. 2006;13:323–30. [PubMed: 16888436]

- Davis KB, Slichter SJ, Corash L. Corrected count increment and percent platelet recovery as measures of posttransfusion platelet response: Problems and a solution. *Transfusion*. 1999;39:586–92. [PubMed: 10378838]

Figures and Tables

Table 1

	Day 1	Day 3	Day 5
Platelet count (10 ¹¹)	3.38	3.36	3.32
pH	7.06	7.03	6.90
Glucose (mg/dl)	353.4	335.1	319.1
Lactate (mmol/L)	11.25	12.69	14.34

Platelet counts and biochemical parameters of BCPP (Average values)

Table 2

	BC pooled platelets	Apheresis platelets
Platelet count (x 10 ¹¹)	3.38	3.30
pH	7.06	6.95

Comparison of BCPP vs Aphaeresis (Average values, Day 1)

Table 3

	Average Platelet increment	**CCI -1 After 2 hour	***PPR % after 2 hour	PPR% at 24 hours
BCPP	32000	7.5	82	46
SDP	34000	7.7	85	52

CCI - Corrected count index, * PPR - Post-transfusion platelet recover

Comparison of efficacy of BCPP and SDP after transfusion

REGISTERED OFFICE

Terumo Penpol Private Limited
P.B.No. 6105, I-2 , Jawahar Nagar,
Thiruvananthapuram - 695003, Kerala, India.
Ph.No : 0471-3015500 / 501
Fax : 0471-2721519
CIN : U33112KL1985PTC004531

www.terumopenpol.com
www.terumobct.com

Toll Free No. 1800 3000 7070

CENTRAL SALES OFFICE

No. 9, Padmanabha View
Gandhi Nagar 2nd Main Road.
Adayar, Chennai - 600 020.
Tel: 044-42054538, 044-30957200
tplmarketing@terumopenpol.com

MARKETING OFFICE

YGRA-12, TC 4/1399(3),
Kuravankonam,
Kowdiar Post,
Thiruvananthapuram-695003
Tel : 0471-3015647/49
pmt@terumopenpol.com

ZONAL OFFICE

DELHI

3E/9, Second Floor,
Jhandewalan Etn, Delhi
Tel : 011-46921062

KOLKATA

839/1, Lake Town, Block A
Kolkatta -700089
West Bengal
Tel : 033-25344609

MUMBAI

Divya Jyot CHS Ltd.,
B-7, Ground floor,
R. K. Singh Marg, Old Nagardas Road,
Andheri (East), Mumbai – 400 069
Tel : 022-2824 0304 / 09320231762

TERUMO PENPOL®