

WB-RP BLOOD BAG SYSTEM

WITH INTEGRAL WHOLE BLOOD FILTER REMOVING PLATELETS

TERUMO PENPOL®

Improve Collections, Streamline Operations

WB-RP Blood Bag System is a unique system designed for collecting and processing whole blood into high-quality, leukodepleted red blood cells and leukodepleted plasma.

■ Reliable Filtration

High performance Sepacell RZ-2000 leukodepletion filter helps create high-quality products that meet local and global standards for room temperature and chilled blood filtration

■ Sharp, dependable phlebotomy needle

Needle with ultra-thin wall and triple bevel design reduces penetration resistance, reduces discomfort of venipuncture and provides stable blood flow

■ Available with comprehensive sampling and safety features

Offers the risk-reducing safety features most important to your blood center


■ Consistent, dependable quality manufacturing

Controlled manufacturing process helps ensure consistent collection, production and storage of high-quality blood components

WB-RP BLOOD BAG SYSTEM


WITH INTEGRAL WHOLE BLOOD FILTER REMOVING PLATELETS

For the production of high-quality, leukodepleted red blood cells and leukodepleted plasma



- Sepacell RZ-2000 Leukodepletion Filter: High performance leukodepletion, high red blood cell recovery, lot traceability
- Diversion Sampling System: Easy to use, air excluded from sample
- Venipuncture Needle: Triple bevel, ultra-thin wall
- Needle Injury Protector: Clear plastic design, simple to use and click confirms needle is locked

WHITE BLOOD CELL FILTRATION PERFORMANCE



Solution Type

Product Code	Collection Volume	Bag Configuration	Description
PB-4W0456J0Y	450 mL	Quadruple	O CPD/SAG-M2

TERUMO PENPOL®

REGISTERED OFFICE
 TERUMO PENPOL Private Limited
 I -2, Jawahar Nagar, Kowdiar P.O,
 Thiruvananthapuram-695003, Kerala, India.
 Ph.No : 0471-3015500 / 501.
 Fax: 0471-2721519
 info@terumopenpol.com

CIN: U33112KL1985PTC004531
www.terumopenpol.com | www.terumobct.com | Toll Free No : 1800 3000 7070

CENTRAL SALES OFFICE
 TERUMO PENPOL Private Limited
 No-9, Padmanabha view – Second Main Road,
 Gandhi Nagar, Adayar, Chennai-600 020
 Tamil Nadu, India.
 Tel: 044-42054538, 42054508
 tplmarketing@terumopenpol.com

MARKETING OFFICE
 TERUMO PENPOL Private Limited
 YGRA-12, TC 4/1399(3),
 Kuravankonam, Kowdiar Post,
 Thiruvananthapuram-695003, Kerala, India.
 Tel: 0471-3015647/649