

Plasma Expressor M1.0



User-Friendly System for the
Manual Separation of Plasma

Plasma Expressor M1.0 Features

- Stainless steel spring-loaded acrylic plate exerts uniform pressure on the blood bag
- Transparent plate allows easy visibility of the bag contents
- Accommodates all types of blood bags up to volumes of 500 mL
- Handle with ball end provides a comfortable grip
- Hook keeps the handle in place before and after expression
- Easy to transport

Plasma Expressor M1.0 Specifications

Mode of operation	Manual
Pressing mechanism	Spring-loaded acrylic plate
Pressing force	3.5 kgf (maximum)
Weight	4.25 kg (approx.)
Finish	Powder coating
Quality standards	ISO 9001
Regulatory compliance	Machinery Directive 2006/42/EC (CE mark)
Product ordering code	M10000M000

Due to continuous development, specifications are subject to change without prior notice. Available in select markets.

TERUMO PENPOL PRIVATE LIMITED
Bhadra Towers, TC 29/2923, Cotton Hill Road,
Vazhuthacaud, Thycaud P.O.,
Thiruvananthapuram, Kerala - 695014, India.
Phone: +91 471-3015500 / 501
E-Mail: Penpol.info@terumobct.com
TERUMOPENPOL.COM
CIN: U33112kL1985PTC004531

